




NFKBIB Antibody

Product Code	CSB-PA919532
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q15653
Immunogen	Fusion protein of Human NFKBIB
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA,WB;ELISA:1:2000-1:10000,WB:1:1000-1:5000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	-20°C, pH7.4 PBS, 0.05% NaN ₃ , 40% Glycerol
Purification Method	Antigen affinity purification
Isotype	IgG
Alias	nuclear factor of kappa light polypeptide gene enhancer in B-cells inhibitor, beta
Immunogen Species	Homo sapiens (Human)
Target Names	NFKBIB

Image



Gel: 10%SDS-PAGE, Lysate: 40 µg, Lane 1-3: Jurkat cells, hela cells, MCF7 cells, Primary antibody: CSB-PA919532(NFKBIB Antibody) at dilution 1/950, Secondary antibody: Goat anti rabbit IgG at 1/8000 dilution, Exposure time: 10 seconds

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.