




# NFKBIB Antibody

<b>Product Code</b>	CSB-PA793544
<b>Storage</b>	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
<b>Uniprot No.</b>	Q15653
<b>Immunogen</b>	Fusion protein of Human NFKBIB
<b>Raised In</b>	Rabbit
<b>Species Reactivity</b>	Human
<b>Tested Applications</b>	ELISA,WB;ELISA:1:2000-1:5000,WB:1:500-1:2000
<b>Form</b>	Liquid
<b>Conjugate</b>	Non-conjugated
<b>Storage Buffer</b>	-20°C, pH7.4 PBS, 0.05% NaN <sub>3</sub> , 40% Glycerol
<b>Purification Method</b>	Antigen affinity purification
<b>Isotype</b>	IgG
<b>Alias</b>	nuclear factor of kappa light polypeptide gene enhancer in B-cells inhibitor, beta
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Target Names</b>	NFKBIB

## Image



Gel: 10%SDS-PAGE, Lysate: 40 µg, Lane 1-3: Jurkat cells, MCF7 cells, hela cells, Primary antibody: CSB-PA793544(NFKBIB Antibody) at dilution 1/1050, Secondary antibody: Goat anti rabbit IgG at 1/8000 dilution, Exposure time: 40 seconds

## Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.