





NFKBID Antibody

Product Code	CSB-PA173796
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8NI38
Immunogen	Fusion protein of Human NFKBID
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA,WB,IHC;ELISA:1:2000-1:5000,WB:1:500-1:2000,IHC:1:25-1:100
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	-20°C, pH7.4 PBS, 0.05% NaN ₃ , 40% Glycerol
Purification Method	Antigen affinity purification
Isotype	IgG
Alias	nuclear factor of kappa light polypeptide gene enhancer in B-cells inhibitor, delta
Immunogen Species	Homo sapiens (Human)
Target Names	NFKBID

Image



The image on the left is immunohistochemistry of paraffin-embedded Human tonsil tissue using CSB-PA173796(NFKBID Antibody) at dilution 1/25, on the right is treated with fusion protein. (Original magnification: x200)



Gel: 10+12%SDS-PAGE, Lysate: 40 µg, Lane 1-2: Hela cells, 231 cells, Primary antibody: CSB-PA173796(NFKBID Antibody) at dilution 1/450, Secondary antibody: Goat anti rabbit IgG at 1/8000 dilution, Exposure time: 2 minutes

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.