




NLRP3 Antibody

| | |
|----------------------------|---|
| Product Code | CSB-PA993403 |
| Storage | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze. |
| Uniprot No. | Q96P20 |
| Immunogen | Fusion protein of Human NLRP3 |
| Raised In | Rabbit |
| Species Reactivity | Human |
| Tested Applications | ELISA,WB;ELISA:1:1000-1:2000,WB:1:200-1:1000 |
| Form | Liquid |
| Conjugate | Non-conjugated |
| Storage Buffer | Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide. |
| Purification Method | Antigen affinity purification |
| Isotype | IgG |
| Product Type | Polyclonal Antibody |
| Immunogen Species | Homo sapiens (Human) |
| Target Names | NLRP3 |

Image



Gel: 8%SDS-PAGE, Lysate: 40 µg, Lane 1-2: K562 cells, hela cells, Primary antibody: CSB-PA993403(NLRP3 Antibody) at dilution 1/88, Secondary antibody: Goat anti rabbit IgG at 1/8000 dilution, Exposure time: 3 minutes

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.