




BAG5 Antibody

Product Code	CSB-PA438344
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9UL15
Immunogen	Fusion protein of Human BAG5
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	ELISA,WB;ELISA:1:500-1:5000,WB:1:500-1:2000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	-20°C, pH7.4 PBS, 0.05% NaN ₃ , 40% Glycerol
Purification Method	Antigen affinity purification
Isotype	IgG
Alias	Bcl-2-associated athanogene 5
Immunogen Species	Homo sapiens (Human)
Target Names	BAG5

Image



Gel: 10%SDS-PAGE, Lysate: 40 µg, Lane 1-3: A549 cells, Hela cells, 293T cells, Primary antibody: CSB-PA438344(BAG5 Antibody) at dilution 1/500, Secondary antibody: Goat anti rabbit IgG at 1/8000 dilution, Exposure time: 1 minute

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.