




# BCOR Antibody

<b>Product Code</b>	CSB-PA186525
<b>Storage</b>	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
<b>Uniprot No.</b>	Q6W2J9
<b>Immunogen</b>	Fusion protein of Human BCOR
<b>Raised In</b>	Rabbit
<b>Species Reactivity</b>	Human, Mouse
<b>Tested Applications</b>	ELISA, WB; ELISA: 1:2000-1:5000, WB: 1:500-1:2000
<b>Form</b>	Liquid
<b>Conjugate</b>	Non-conjugated
<b>Storage Buffer</b>	-20°C, pH7.4 PBS, 0.05% NaN <sub>3</sub> , 40% Glycerol
<b>Purification Method</b>	Antigen affinity purification
<b>Isotype</b>	IgG
<b>Alias</b>	BCL6 corepressor
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Target Names</b>	BCOR

## Image



Gel: 8%SDS-PAGE, Lysate: 40 µg, Lane: Mouse stomach tissue, Primary antibody: CSB-PA186525(BCOR Antibody) at dilution 1/350, Secondary antibody: Goat anti rabbit IgG at 1/8000 dilution, Exposure time: 40 seconds

## Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.