





BCL2L12 Antibody

Product Code	CSB-PA582788
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9HB09
Immunogen	Fusion protein of Human BCL2L12
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA,IHC;ELISA:1:1000-1:2000,IHC:1:25-1:100
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	-20°C, pH7.4 PBS, 0.05% NaN ₃ , 40% Glycerol
Purification Method	Antigen affinity purification
Isotype	IgG
Alias	BCL2-like 12 (proline rich)
Immunogen Species	Homo sapiens (Human)
Target Names	BCL2L12

Image



The image is immunohistochemistry of paraffin-embedded Human thyroid cancer tissue using CSB-PA582788(BCL2L12 Antibody) at dilution 1/20. (Original magnification: x200)



The image is immunohistochemistry of paraffin-embedded Human gastric cancer tissue using CSB-PA582788(BCL2L12 Antibody) at dilution 1/20. (Original magnification: x200)

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.