





BCL2A1 Antibody

Product Code	CSB-PA561279
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q16548
Immunogen	Fusion protein of Human BCL2A1
Raised In	Rabbit
Species Reactivity	Human,Mouse
Tested Applications	ELISA,WB,IHC;ELISA:1:2000-1:5000,WB:1:200-1:1000,IHC:1:25-1:100
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	-20°C, pH7.4 PBS, 0.05% NaN ₃ , 40% Glycerol
Purification Method	Antigen affinity purification
Isotype	IgG
Alias	BCL2-related protein A1
Immunogen Species	Homo sapiens (Human)
Target Names	BCL2A1

Image



The image on the left is immunohistochemistry of paraffin-embedded Human thyroid cancer tissue using CSB-PA561279(BCL2A1 Antibody) at dilution 1/20, on the right is treated with fusion protein. (Original magnification: x200)



Gel: 12%SDS-PAGE, Lysate: 40 µg, Lane 1-2: Mouse testis tissue, hela cells, Primary antibody: CSB-PA561279(BCL2A1 Antibody) at dilution 1/117, Secondary antibody: Goat anti rabbit IgG at 1/8000 dilution, Exposure time: 1 minute