





IL10RB Antibody

Product Code	CSB-PA294228
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q08334
Immunogen	Fusion protein of Human IL10RB
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA,IHC;ELISA:1:1000-1:10000,IHC:1:25-1:100
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	-20°C, pH7.4 PBS, 0.05% NaN ₃ , 40% Glycerol
Purification Method	Antigen affinity purification
Isotype	IgG
Alias	Interleukin-10 receptor subunit beta
Immunogen Species	Homo sapiens (Human)
Target Names	IL10RB

Image



The image is immunohistochemistry of paraffin-embedded Human cervical cancer tissue using CSB-PA294228(IL10RB Antibody) at dilution 1/20. (Original magnification: ×200)



The image is immunohistochemistry of paraffin-embedded Human esophagus cancer tissue using CSB-PA294228(IL10RB Antibody) at dilution 1/20. (Original magnification: ×200)