





IL13RA2 Antibody

Product Code	CSB-PA239289
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q14627
Immunogen	Fusion protein of Human IL13RA2
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA,IHC;ELISA:1:1000-1:2000,IHC:1:15-1:50
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	-20°C, pH7.4 PBS, 0.05% NaN ₃ , 40% Glycerol
Purification Method	Antigen affinity purification
Isotype	IgG
Alias	interleukin 13 receptor, alpha 2
Immunogen Species	Homo sapiens (Human)
Target Names	IL13RA2

Image



The image on the left is immunohistochemistry of paraffin-embedded Human lymphoma tissue using CSB-PA239289(IL13RA2 Antibody) at dilution 1/30, on the right is treated with fusion protein. (Original magnification: x200)



The image on the left is immunohistochemistry of paraffin-embedded Human cervical cancer tissue using CSB-PA239289(IL13RA2 Antibody) at dilution 1/30, on the right is treated with fusion protein. (Original magnification: x200)

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.