




IL17RD Antibody

| | |
|----------------------------|---|
| Product Code | CSB-PA206405 |
| Storage | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze. |
| Uniprot No. | Q8NFM7 |
| Immunogen | Fusion protein of Human IL17RD |
| Raised In | Rabbit |
| Species Reactivity | Human, Mouse |
| Tested Applications | ELISA, IHC; ELISA: 1:500-1:2000, IHC: 1:5-1:20 |
| Form | Liquid |
| Conjugate | Non-conjugated |
| Storage Buffer | -20°C, pH7.4 PBS, 0.05% NaN ₃ , 40% Glycerol |
| Purification Method | Antigen affinity purification |
| Isotype | IgG |
| Alias | Interleukin-17 receptor D |
| Immunogen Species | Homo sapiens (Human) |
| Target Names | IL17RD |

Image



The image on the left is immunohistochemistry of paraffin-embedded Human lung cancer tissue using CSB-PA206405 (IL17RD Antibody) at dilution 1/5, on the right is treated with fusion protein. (Original magnification: x200)

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.