




IL1R2 Antibody

Product Code	CSB-PA247988
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P27930
Immunogen	Fusion protein of Human IL1R2
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA,IHC;ELISA:1:1000-1:5000,IHC:1:25-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	-20°C, pH7.4 PBS, 0.05% NaN ₃ , 40% Glycerol
Purification Method	Antigen affinity purification
Isotype	IgG
Alias	Interleukin-1 receptor type 2
Immunogen Species	Homo sapiens (Human)
Target Names	IL1R2

Image



The image on the left is immunohistochemistry of paraffin-embedded Human lung cancer tissue using CSB-PA247988(IL1R2 Antibody) at dilution 1/20, on the right is treated with fusion protein. (Original magnification: x200)