





IL20RA Antibody

Product Code	CSB-PA922324
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9UHF4
Immunogen	Fusion protein of Human IL20RA
Raised In	Rabbit
Species Reactivity	Human, Mouse
Tested Applications	ELISA, IHC; ELISA: 1:2000-1:5000, IHC: 1:50-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	-20°C, pH7.4 PBS, 0.05% NaN ₃ , 40% Glycerol
Purification Method	Antigen affinity purification
Isotype	IgG
Alias	interleukin 20 receptor, alpha
Immunogen Species	Homo sapiens (Human)
Target Names	IL20RA

Image



The image on the left is immunohistochemistry of paraffin-embedded Human colon cancer tissue using CSB-PA922324(IL20RA Antibody) at dilution 1/35, on the right is treated with fusion protein. (Original magnification: x200)



The image on the left is immunohistochemistry of paraffin-embedded Human lung cancer tissue using CSB-PA922324(IL20RA Antibody) at dilution 1/35, on the right is treated with fusion protein. (Original magnification: x200)

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.