





IL23R Antibody

Product Code	CSB-PA273292
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q5VWK5
Immunogen	Fusion protein of Human IL23R
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, WB, IHC; ELISA: 1:500-1:5000, WB: 1:500-1:2000, IHC: 1:5-1:20
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	-20°C, pH7.4 PBS, 0.05% NaN ₃ , 40% Glycerol
Purification Method	Antigen affinity purification
Isotype	IgG
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	IL23R

Image



The image on the left is immunohistochemistry of paraffin-embedded Human rectal cancer tissue using CSB-PA273292(IL23R Antibody) at dilution 1/10, on the right is treated with fusion protein. (Original magnification: x200)



Gel: 8%SDS-PAGE, Lysate: 40 µg, Lane: Human sigmoid tissue, Primary antibody: CSB-PA273292(IL23R Antibody) at dilution 1/300, Secondary antibody: Goat anti rabbit IgG at 1/8000 dilution, Exposure time: 30 seconds

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.