





HTR2B Antibody

Product Code	CSB-PA092212
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P41595
Immunogen	Fusion protein of Human HTR2B
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, WB, IHC; ELISA: 1:500-1:2000, WB: 1:200-1:1000, IHC: 1:25-1:100
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	-20°C, pH7.4 PBS, 0.05% NaN ₃ , 40% Glycerol
Purification Method	Antigen affinity purification
Isotype	IgG
Alias	5-hydroxytryptamine receptor 2B
Immunogen Species	Homo sapiens (Human)
Target Names	HTR2B

Image



The image on the left is immunohistochemistry of paraffin-embedded Human lung cancer tissue using CSB-PA092212(HTR2B Antibody) at dilution 1/20, on the right is treated with fusion protein. (Original magnification: x200)



Gel: 10%SDS-PAGE, Lysate: 20 µg, Lane: 293T cells, Primary antibody: CSB-PA092212(HTR2B Antibody) at dilution 1/500, Secondary antibody: Goat anti rabbit IgG at 1/8000 dilution, Exposure time: 1 minute