




ABCB5 Antibody

Product Code	CSB-PA967715
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q2M3G0
Immunogen	Fusion protein of Human ABCB5
Raised In	Rabbit
Species Reactivity	Human,Mouse
Tested Applications	ELISA,IHC;ELISA:1:500-1:2000,IHC:1:10-1:50
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	-20°C, pH7.4 PBS, 0.05% NaN ₃ , 40% Glycerol
Purification Method	Antigen affinity purification
Isotype	IgG
Alias	ATP-binding cassette sub-family B member 5
Immunogen Species	Homo sapiens (Human)
Target Names	ABCB5

Image



The image on the left is immunohistochemistry of paraffin-embedded Human lung cancer tissue using CSB-PA967715(ABCB5 Antibody) at dilution 1/10, on the right is treated with fusion protein. (Original magnification: x200)