





ANGPTL4 Antibody

Product Code	CSB-PA261295
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9BY76
Immunogen	Fusion protein of Human ANGPTL4
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA,IHC;ELISA:1:2000-1:5000,IHC:1:25-1:100
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	-20°C, pH7.4 PBS, 0.05% NaN ₃ , 40% Glycerol
Purification Method	Antigen affinity purification
Isotype	IgG
Alias	Angiopietin-related protein 4
Immunogen Species	Homo sapiens (Human)
Target Names	ANGPTL4

Image



The image is immunohistochemistry of paraffin-embedded Human cervical cancer tissue using CSB-PA261295(ANGPTL4 Antibody) at dilution 1/25. (Original magnification: ×200)



The image is immunohistochemistry of paraffin-embedded Human esophagus cancer tissue using CSB-PA261295(ANGPTL4 Antibody) at dilution 1/25. (Original magnification: ×200)

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.