





ANXA8 Antibody

Product Code	CSB-PA944795
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P13928
Immunogen	Fusion protein of Human ANXA8
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	ELISA,IHC;ELISA:1:2000-1:5000,IHC:1:25-1:100
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	-20°C, pH7.4 PBS, 0.05% NaN ₃ , 40% Glycerol
Purification Method	Antigen affinity purification
Isotype	IgG
Alias	annexin A8
Immunogen Species	Homo sapiens (Human)
Target Names	ANXA8

Image



The image on the left is immunohistochemistry of paraffin-embedded Human brain tissue using CSB-PA944795(ANXA8 Antibody) at dilution 1/25, on the right is treated with fusion protein. (Original magnification: ×200)



The image on the left is immunohistochemistry of paraffin-embedded Human lung cancer tissue using CSB-PA944795(ANXA8 Antibody) at dilution 1/25, on the right is treated with fusion protein. (Original magnification: ×200)