





KRT23 Antibody

Product Code	CSB-PA041845
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9C075
Immunogen	Fusion protein of Human KRT23
Raised In	Rabbit
Species Reactivity	Human, Mouse
Tested Applications	ELISA, WB, IHC; ELISA: 1:2000-1:5000, WB: 1:500-1:2000, IHC: 1:25-1:100
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	-20°C, pH7.4 PBS, 0.05% NaN ₃ , 40% Glycerol
Purification Method	Antigen affinity purification
Isotype	IgG
Alias	keratin 23 (histone deacetylase inducible)
Immunogen Species	Homo sapiens (Human)
Target Names	KRT23

Image



The image on the left is immunohistochemistry of paraffin-embedded Human colon cancer tissue using CSB-PA041845(KRT23 Antibody) at dilution 1/20, on the right is treated with fusion protein. (Original magnification: x200)



Gel: 10%SDS-PAGE, Lysate: 30 μg, Lane: HeLa cells, Primary antibody: CSB-PA041845(KRT23 Antibody) at dilution 1/300, Secondary antibody: Goat anti rabbit IgG at 1/8000 dilution, Exposure time: 20 seconds

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.