





TTR Antibody

Product Code	CSB-PA143118
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P02766
Immunogen	Fusion protein of Human TTR
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	ELISA,IHC;ELISA:1:1000-1:5000,IHC:1:25-1:100
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	-20°C, pH7.4 PBS, 0.05% NaN ₃ , 40% Glycerol
Purification Method	Antigen affinity purification
Isotype	IgG
Alias	Transthyretin
Immunogen Species	Homo sapiens (Human)
Target Names	TTR

Image



The image on the left is immunohistochemistry of paraffin-embedded Human stomach cancer tissue using CSB-PA143118(TTR Antibody) at dilution 1/10, on the right is treated with fusion protein. (Original magnification: x200)



The image on the left is immunohistochemistry of paraffin-embedded Human tonsil tissue using CSB-PA143118(TTR Antibody) at dilution 1/10, on the right is treated with fusion protein. (Original magnification: x200)

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.