




# BSA Antibody

<b>Product Code</b>	CSB-PA592874
<b>Storage</b>	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
<b>Uniprot No.</b>	P02769
<b>Immunogen</b>	Fusion protein of BSA
<b>Raised In</b>	Rabbit
<b>Species Reactivity</b>	Bovine
<b>Tested Applications</b>	ELISA,WB;ELISA:1:2000-1:10000,WB:1:1000-1:5000
<b>Form</b>	Liquid
<b>Conjugate</b>	Non-conjugated
<b>Storage Buffer</b>	-20°C, pH7.4 PBS, 0.05% NaN <sub>3</sub> , 40% Glycerol
<b>Purification Method</b>	Antigen affinity purification
<b>Isotype</b>	IgG
<b>Alias</b>	Bovine serum albumin
<b>Immunogen Species</b>	Bos taurus (Bovine)
<b>Target Names</b>	ALB

## Image



Gel: 8%SDS-PAGE, Lysate: 0.05-0.4 µg, Lane 1-4: 0.05 µg, 0.1 µg, 0.2 µg, 0.4 µg BSA, Primary antibody: CSB-PA592874(BSA Antibody) at dilution 1/1450, Secondary antibody: Goat anti rabbit IgG at 1/8000 dilution, Exposure time: 1 second

## Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.