





# ATP5J Antibody

<b>Product Code</b>	CSB-PA075471
<b>Storage</b>	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
<b>Uniprot No.</b>	P18859
<b>Immunogen</b>	Fusion protein of Human ATP5J
<b>Raised In</b>	Rabbit
<b>Species Reactivity</b>	Human, Mouse
<b>Tested Applications</b>	ELISA, WB, IHC; ELISA: 1:1000-1:2000, WB: 1:200-1:1000, IHC: 1:15-1:50
<b>Form</b>	Liquid
<b>Conjugate</b>	Non-conjugated
<b>Storage Buffer</b>	-20°C, pH7.4 PBS, 0.05% NaN <sub>3</sub> , 40% Glycerol
<b>Purification Method</b>	Antigen affinity purification
<b>Isotype</b>	IgG
<b>Alias</b>	ATP synthase, H <sup>+</sup> transporting, mitochondrial Fo complex, subunit F6
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Target Names</b>	ATP5J


## Image



The image on the left is immunohistochemistry of paraffin-embedded Human brain tissue using CSB-PA075471(ATP5J Antibody) at dilution 1/15, on the right is treated with fusion protein. (Original magnification: x200)



The image on the left is immunohistochemistry of paraffin-embedded Human thyroid cancer tissue using CSB-PA075471(ATP5J Antibody) at dilution 1/15, on the right is treated with fusion protein. (Original magnification: x200)



Gel: 12%SDS-PAGE, Lysate: 40  $\mu$ g, Lane: Lovo cells, Primary antibody: CSB-PA075471(ATP5J Antibody) at dilution 1/200, Secondary antibody: Goat anti rabbit IgG at 1/8000 dilution, Exposure time: 5 minutes

## Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.