




BCL2L1 Antibody

Product Code	CSB-PA223208
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q07817
Immunogen	Fusion protein of Human BCL2L1
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	ELISA,WB;ELISA:1:2000-1:5000,WB:1:500-1:2000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	-20°C, pH7.4 PBS, 0.05% NaN ₃ , 40% Glycerol
Purification Method	Antigen affinity purification
Isotype	IgG
Alias	Bcl-2-like protein 1
Immunogen Species	Homo sapiens (Human)
Target Names	BCL2L1

Image



Gel: 12+15%SDS-PAGE, Lysate: 30 µg, Lane: Human cervical cancer tissue, Primary antibody: CSB-PA223208(BCL2L1 Antibody) at dilution 1/450, Secondary antibody: Goat anti rabbit IgG at 1/8000 dilution, Exposure time: 10 seconds

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.