





TNFSF13B Antibody

Product Code	CSB-PA108957
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9Y275
Immunogen	Fusion protein of Human TNFSF13B
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA,IHC;ELISA:1:1000-1:5000,IHC:1:50-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	-20°C, pH7.4 PBS, 0.05% NaN ₃ , 40% Glycerol
Purification Method	Antigen affinity purification
Isotype	IgG
Alias	tumor necrosis factor (ligand) superfamily, member 13b
Immunogen Species	Homo sapiens (Human)
Target Names	TNFSF13B

Image



The image on the left is immunohistochemistry of paraffin-embedded Human liver cancer tissue using CSB-PA108957(TNFSF13B Antibody) at dilution 1/30, on the right is treated with fusion protein. (Original magnification: x200)



The image on the left is immunohistochemistry of paraffin-embedded Human lymphoma tissue using CSB-PA108957(TNFSF13B Antibody) at dilution 1/30, on the right is treated with fusion protein. (Original magnification: x200)

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.