





BID Antibody

Product Code	CSB-PA787438
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P55957
Immunogen	Fusion protein of Human BID
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA,WB,IHC;ELISA:1:1000-1:5000,WB:1:200-1:2000,IHC:1:25-1:100
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	-20°C, pH7.4 PBS, 0.05% NaN ₃ , 40% Glycerol
Purification Method	Antigen affinity purification
Isotype	IgG
Alias	BH3-interacting domain death agonist
Immunogen Species	Homo sapiens (Human)
Target Names	BID

Image



The image on the left is immunohistochemistry of paraffin-embedded Human renal cancer tissue using CSB-PA787438(BID Antibody) at dilution 1/25, on the right is treated with fusion protein. (Original magnification: x200)



Gel: 12%SDS-PAGE, Lysate: 30 µg, Lane 1-4: Human fetal brain tissue, Human liver tissue, HT29 cells, HepG2 cells, Primary antibody: CSB-PA787438(BID Antibody) at dilution 1/260, Secondary antibody: Goat anti rabbit IgG at 1/8000 dilution, Exposure time: 10 seconds

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.