




# CAPN2 Antibody

<b>Product Code</b>	CSB-PA257142
<b>Storage</b>	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
<b>Uniprot No.</b>	P17655
<b>Immunogen</b>	Fusion protein of Human CAPN2
<b>Raised In</b>	Rabbit
<b>Species Reactivity</b>	Human,Mouse,Rat
<b>Tested Applications</b>	ELISA,WB;ELISA:1:1000-1:5000,WB:1:200-1:1000
<b>Form</b>	Liquid
<b>Conjugate</b>	Non-conjugated
<b>Storage Buffer</b>	-20°C, pH7.4 PBS, 0.05% NaN <sub>3</sub> , 40% Glycerol
<b>Purification Method</b>	Antigen affinity purification
<b>Isotype</b>	IgG
<b>Alias</b>	Calpain-2 catalytic subunit
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Target Names</b>	CAPN2

## Image



Gel: 10%SDS-PAGE, Lysate: 30 µg, Lane 1-3: Mouse brain tissue, Mouse skin tissue, 293T cells, Primary antibody: CSB-PA257142(CAPN2 Antibody) at dilution 1/200, Secondary antibody: Goat anti rabbit IgG at 1/8000 dilution, Exposure time: 40 seconds

## Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.