




CAPN2 Antibody

Product Code	CSB-PA257142
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P17655
Immunogen	Fusion protein of Human CAPN2
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	ELISA,WB;ELISA:1:1000-1:5000,WB:1:200-1:1000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	-20°C, pH7.4 PBS, 0.05% NaN ₃ , 40% Glycerol
Purification Method	Antigen affinity purification
Isotype	IgG
Alias	Calpain-2 catalytic subunit
Immunogen Species	Homo sapiens (Human)
Target Names	CAPN2

Image



Gel: 10%SDS-PAGE, Lysate: 30 µg, Lane 1-3: Mouse brain tissue, Mouse skin tissue, 293T cells, Primary antibody: CSB-PA257142(CAPN2 Antibody) at dilution 1/200, Secondary antibody: Goat anti rabbit IgG at 1/8000 dilution, Exposure time: 40 seconds

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.