




CA1 Antibody

Product Code	CSB-PA348998
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P00915
Immunogen	Fusion protein of Human CA1
Raised In	Rabbit
Species Reactivity	Human,Rat
Tested Applications	ELISA,WB;ELISA:1:1000-1:5000,WB:1:200-1:1000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	-20°C, pH7.4 PBS, 0.05% NaN ₃ , 40% Glycerol
Purification Method	Antigen affinity purification
Isotype	IgG
Alias	Carbonic anhydrase 1
Immunogen Species	Homo sapiens (Human)
Target Names	CA1

Image



Gel: 10%SDS-PAGE, Lysate: 24 µg, Lane:
Human lung cancer tissue, Primary antibody:
CSB-PA348998(CA1 Antibody) at dilution 1/350,
Secondary antibody: Goat anti rabbit IgG at
1/8000 dilution, Exposure time: 1 minute

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.