




CA14 Antibody

Product Code	CSB-PA700069
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9ULX7
Immunogen	Fusion protein of Human CA14
Raised In	Rabbit
Species Reactivity	Human,Mouse
Tested Applications	ELISA,WB;ELISA:1:2000-1:10000,WB:1:1000-1:5000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	-20°C, pH7.4 PBS, 0.05% NaN ₃ , 40% Glycerol
Purification Method	Antigen affinity purification
Isotype	IgG
Alias	Carbonic anhydrase XIV
Immunogen Species	Homo sapiens (Human)
Target Names	CA14

Image



Gel: 10%SDS-PAGE, Lysate: 30 µg, Lane 1-2: Human liver cancer tissue, mouse heart tissue, Primary antibody: CSB-PA700069(CA14 Antibody) at dilution 1/1050, Secondary antibody: Goat anti rabbit IgG at 1/8000 dilution, Exposure time: 20 seconds

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.