





CTNNA1 Antibody

Product Code	CSB-PA218040
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P35221
Immunogen	Fusion protein of Human CTNNA1
Raised In	Rabbit
Species Reactivity	Human, Mouse
Tested Applications	ELISA, IHC; ELISA: 1:2000-1:5000, IHC: 1:25-1:100
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	-20°C, pH7.4 PBS, 0.05% NaN ₃ , 40% Glycerol
Purification Method	Antigen affinity purification
Isotype	IgG
Alias	Catenin (cadherin-associated protein), alpha 1, 102kDa
Immunogen Species	Homo sapiens (Human)
Target Names	CTNNA1

Image



The image on the left is immunohistochemistry of paraffin-embedded Human brain tissue using CSB-PA218040(CTNNA1 Antibody) at dilution 1/15, on the right is treated with fusion protein. (Original magnification: x200)



The image on the left is immunohistochemistry of paraffin-embedded Human thyroid cancer tissue using CSB-PA218040(CTNNA1 Antibody) at dilution 1/15, on the right is treated with fusion protein. (Original magnification: x200)

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.