





CTSL Antibody

Product Code	CSB-PA092272
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P07711
Immunogen	Fusion protein of Human CTSL
Raised In	Rabbit
Species Reactivity	Human,Mouse
Tested Applications	ELISA,WB,IHC;ELISA:1:2000-1:5000,WB:1:500-1:2000,IHC:1:50-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	-20°C, pH7.4 PBS, 0.05% NaN ₃ , 40% Glycerol
Purification Method	Antigen affinity purification
Isotype	IgG
Alias	cathepsin L
Immunogen Species	Homo sapiens (Human)
Target Names	CTSL


Image



The image on the left is immunohistochemistry of paraffin-embedded Human brain tissue using CSB-PA092272(CTSL Antibody) at dilution 1/40, on the right is treated with fusion protein. (Original magnification: x200)



The image on the left is immunohistochemistry of paraffin-embedded Human thyroid cancer tissue using CSB-PA092272(CTSL Antibody) at dilution 1/40, on the right is treated with fusion protein. (Original magnification: x200)



Gel: 15%SDS-PAGE, Lysate: 40 μ g, Lane 1-2:
 Mouse kidney tissue, A172 cells, Primary
 antibody: CSB-PA092272(CTSL Antibody) at
 dilution 1/900, Secondary antibody: Goat anti
 rabbit IgG at 1/8000 dilution, Exposure time: 10
 seconds