





CBX8 Antibody

Product Code	CSB-PA557525
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9HC52
Immunogen	Fusion protein of Human CBX8
Raised In	Rabbit
Species Reactivity	Human, Mouse
Tested Applications	ELISA, WB, IHC; ELISA: 1:2000-1:5000, WB: 1:500-1:2000, IHC: 1:10-1:50
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	-20°C, pH7.4 PBS, 0.05% NaN ₃ , 40% Glycerol
Purification Method	Antigen affinity purification
Isotype	IgG
Alias	chromobox homolog 8
Immunogen Species	Homo sapiens (Human)
Target Names	CBX8

Image



The image on the left is immunohistochemistry of paraffin-embedded Human lung cancer tissue using CSB-PA557525(CBX8 Antibody) at dilution 1/10, on the right is treated with fusion protein. (Original magnification: x200)



Gel: 10%SDS-PAGE, Lysate: 40 µg, Lane 1-3: 231 cells, HeLa cells, A549 cells, Primary antibody: CSB-PA557525(CBX8 Antibody) at dilution 1/200, Secondary antibody: Goat anti rabbit IgG at 1/8000 dilution, Exposure time: 1 minute

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.