




CBX5 Antibody

Product Code	CSB-PA914989
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P45973
Immunogen	Fusion protein of Human CBX5
Raised In	Rabbit
Species Reactivity	Human,Mouse
Tested Applications	ELISA,WB;ELISA:1:1000-1:5000,WB:1:200-1:1000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	-20°C, pH7.4 PBS, 0.05% NaN ₃ , 40% Glycerol
Purification Method	Antigen affinity purification
Isotype	IgG
Alias	Chromobox protein homolog 5
Immunogen Species	Homo sapiens (Human)
Target Names	CBX5

Image



Gel: 10%SDS-PAGE, Lysate: 24 µg, Lane 1-2: HeLa cells, Human endometrial cancer tissue, Primary antibody: CSB-PA914989(CBX5 Antibody) at dilution 1/300, Secondary antibody: Goat anti rabbit IgG at 1/8000 dilution, Exposure time: 1 minute

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.