




CCL13 Antibody

Product Code	CSB-PA218529
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q99616
Immunogen	Fusion protein of Human CCL13
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA,WB;ELISA:1:2000-1:5000,WB:1:500-1:2000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	-20°C, pH7.4 PBS, 0.05% NaN ₃ , 40% Glycerol
Purification Method	Antigen affinity purification
Isotype	IgG
Alias	chemokine (C-C motif) ligand 13
Immunogen Species	Homo sapiens (Human)
Target Names	CCL13

Image



Gel: 12%SDS-PAGE, Lysate: 50 µg, Lane: 293T cells, Primary antibody: CSB-PA218529(CCL13 Antibody) at dilution 1/900, Secondary antibody: Goat anti rabbit IgG at 1/8000 dilution, Exposure time: 10 seconds

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.