





CDK5 Antibody

Product Code	CSB-PA280305
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q00535
Immunogen	Fusion protein of Human CDK5
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	ELISA,WB,IHC;ELISA:1:1000-1:2000,WB:1:200-1:1000,IHC:1:10-1:50
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	-20°C, pH7.4 PBS, 0.05% NaN ₃ , 40% Glycerol
Purification Method	Antigen affinity purification
Isotype	IgG
Alias	cyclin-dependent kinase 5
Immunogen Species	Homo sapiens (Human)
Target Names	CDK5

Image



The image on the left is immunohistochemistry of paraffin-embedded Human liver cancer tissue using CSB-PA280305(CDK5 Antibody) at dilution 1/10, on the right is treated with fusion protein. (Original magnification: x200)



Gel: 10+12%SDS-PAGE, Lysate: 30 µg, Lane 1-2: Human fetal brain tissue, Human colon cancer tissue, Primary antibody: CSB-PA280305(CDK5 Antibody) at dilution 1/150, Secondary antibody: Goat anti rabbit IgG at 1/8000 dilution, Exposure time: 5 minutes

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.