





CDK7 Antibody

Product Code	CSB-PA440881
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P50613
Immunogen	Fusion protein of Human CDK7
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	ELISA,WB,IHC;ELISA:1:1000-1:5000,WB:1:500-1:2000,IHC:1:15-1:50
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	-20°C, pH7.4 PBS, 0.05% NaN ₃ , 40% Glycerol
Purification Method	Antigen affinity purification
Isotype	IgG
Alias	cyclin-dependent kinase 7
Immunogen Species	Homo sapiens (Human)
Target Names	CDK7

Image



The image on the left is immunohistochemistry of paraffin-embedded Human lung cancer tissue using CSB-PA440881(CDK7 Antibody) at dilution 1/30, on the right is treated with fusion protein. (Original magnification: X200)



Gel: 10%SDS-PAGE, Lysate: 60μg, Lane: A549 cells, Primary antibody: CSB-PA440881(CDK7 Antibody) at dilution 1/500, Secondary antibody: Goat anti rabbit IgG at 1/8000 dilution, Exposure time: 2 minutes