





CDKN1A Antibody

Product Code	CSB-PA934176
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P38936
Immunogen	Fusion protein of Human CDKN1A
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, WB, IHC; ELISA: 1:1000-1:5000, WB: 1:500-1:2000, IHC: 1:15-1:50
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	-20°C, pH7.4 PBS, 0.05% NaN ₃ , 40% Glycerol
Purification Method	Antigen affinity purification
Isotype	IgG
Alias	cyclin-dependent kinase inhibitor 1A (p21, Cip1)
Immunogen Species	Homo sapiens (Human)
Target Names	CDKN1A

Image



The image on the left is immunohistochemistry of paraffin-embedded Human lung cancer tissue using CSB-PA934176(CDKN1A Antibody) at dilution 1/15, on the right is treated with fusion protein. (Original magnification: x200)



Gel: 10%SDS-PAGE, Lysate: 40 µg, Lane 1-3: 293T cells, human brain malignant glioma tissue, Human liver cancer tissue, Primary antibody: CSB-PA934176(CDKN1A Antibody) at dilution 1/250, Secondary antibody: Goat anti rabbit IgG at 1/8000 dilution, Exposure time: 5 minutes

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.