





EMC8 Antibody

| | |
|----------------------------|---|
| Product Code | CSB-PA939763 |
| Storage | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze. |
| Uniprot No. | O43402 |
| Immunogen | Full length fusion protein |
| Raised In | Rabbit |
| Species Reactivity | Human,Mouse,Rat |
| Tested Applications | ELISA,WB,IHC;ELISA:1:2000-1:5000,WB:1:500-1:2000,IHC:1:50-1:200 |
| Form | Liquid |
| Conjugate | Non-conjugated |
| Storage Buffer | -20°C, pH7.4 PBS, 0.05% NaN ₃ , 40% Glycerol |
| Purification Method | Antigen affinity purification |
| Isotype | IgG |
| Alias | ER membrane protein complex subunit 8 |
| Immunogen Species | Homo sapiens (Human) |
| Target Names | EMC8 |


Image



The image on the left is immunohistochemistry of paraffin-embedded Human cervical cancer tissue using CSB-PA939763(EMC8 Antibody) at dilution 1/30, on the right is treated with fusion protein. (Original magnification: x200)



The image on the left is immunohistochemistry of paraffin-embedded Human brain tissue using CSB-PA939763(EMC8 Antibody) at dilution 1/30, on the right is treated with fusion protein. (Original magnification: x200)



Gel: 10%SDS-PAGE,Lysate: 40 μ g,Lane 1-2:
Hela cells, HT-29 cells,Primary antibody: CSB-
PA939763(EMC8 Antibody) at dilution 1/650
dilution,Secondary antibody: Goat anti rabbit IgG
at 1/8000 dilution,Exposure time: 1 minute

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.