




CREB1 Antibody

Product Code	CSB-PA205203
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P16220
Immunogen	Fusion protein of Human CREB1
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	ELISA,WB;ELISA:1:1000-1:2000,WB:1:200-1:1000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	-20°C, pH7.4 PBS, 0.05% NaN ₃ , 40% Glycerol
Purification Method	Antigen affinity purification
Isotype	IgG
Alias	cAMP responsive element binding protein 1
Immunogen Species	Homo sapiens (Human)
Target Names	CREB1

Image



Gel: 8%SDS-PAGE, Lysate: 40 µg, Lane:
 NIH/3T3 cells, Primary antibody: CSB-PA205203(CREB1 Antibody) at dilution 1/350,
 Secondary antibody: Goat anti rabbit IgG at 1/8000 dilution, Exposure time: 20 seconds

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.