




CD151 Antibody

Product Code	CSB-PA940948
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P48509
Immunogen	Fusion protein of Human CD151
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	ELISA,IHC;ELISA:1:2000-1:10000,IHC:1:50-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	-20°C, pH7.4 PBS, 0.05% NaN ₃ , 40% Glycerol
Purification Method	Antigen affinity purification
Isotype	IgG
Alias	Cluster of Differentiation 151 antigen
Immunogen Species	Homo sapiens (Human)
Target Names	CD151

Image



The image on the left is immunohistochemistry of paraffin-embedded Human liver cancer tissue using CSB-PA940948(CD151 Antibody) at dilution 1/25, on the right is treated with fusion protein. (Original magnification: x200)

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.