




CUL7 Antibody

Product Code	CSB-PA051064
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q14999
Immunogen	Fusion protein of Human CUL7
Raised In	Rabbit
Species Reactivity	Human,Rat
Tested Applications	ELISA,IHC,WB;ELISA:1:1000-1:2000,IHC:1:25-1:100,WB:1:200-1:1000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	-20°C, pH7.4 PBS, 0.05% NaN ₃ , 40% Glycerol
Purification Method	Antigen affinity purification
Isotype	IgG
Alias	cullin 7
Immunogen Species	Homo sapiens (Human)
Target Names	CUL7

Image



The image is immunohistochemistry of paraffin-embedded human thyroid cancer tissue using CSB-PA051064(CUL7 Antibody) at dilution 1/35. (Original magnification: x200)

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.