




IL18R1 Antibody

Product Code	CSB-PA056452
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q13478
Immunogen	Fusion protein of Human IL18R1
Raised In	Rabbit
Species Reactivity	Human,Mouse
Tested Applications	ELISA,WB;ELISA:1:2000-1:5000,WB:1:500-1:2000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	-20°C, pH7.4 PBS, 0.05% NaN ₃ , 40% Glycerol
Purification Method	Antigen affinity purification
Isotype	IgG
Alias	interleukin 18 receptor 1
Immunogen Species	Homo sapiens (Human)
Target Names	IL18R1

Image



Gel: 8%SDS-PAGE,Lysate: 40 µg,,Primary antibody: CSB-PA056452(IL18R1 Antibody) at dilution 1/200 dilution,Secondary antibody: Goat anti rabbit IgG at 1/8000 dilution,Exposure time: 1 minute