




CYP1A1 Antibody

Product Code	CSB-PA231355
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P04798
Immunogen	Fusion protein of Human CYP1A1
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	ELISA,WB;ELISA:1:1000-1:5000,WB:1:500-1:2000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	-20°C, pH7.4 PBS, 0.05% NaN ₃ , 40% Glycerol
Purification Method	Antigen affinity purification
Isotype	IgG
Alias	cytochrome P450, family 1, subfamily A, polypeptide 1
Immunogen Species	Homo sapiens (Human)
Target Names	CYP1A1

Image



Gel: 10%SDS-PAGE, Lysate: 30 µg, Lane: Mouse liver tissue, Primary antibody: CSB-PA231355(CYP1A1 Antibody) at dilution 1/500, Secondary antibody: Goat anti rabbit IgG at 1/8000 dilution, Exposure time: 5 seconds

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.