




TNFRSF12A Antibody

Product Code	CSB-PA199353
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9NP84
Immunogen	Fusion protein of Human TNFRSF12A
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA,IHC; ELISA:1:1000-1:2000,IHC:1:5-1:20
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	-20°C, pH7.4 PBS, 0.05% NaN ₃ , 40% Glycerol
Purification Method	Antigen affinity purification
Isotype	IgG
Alias	Tumor necrosis factor receptor superfamily, member 12A
Immunogen Species	Homo sapiens (Human)
Target Names	TNFRSF12A

Image



The image on the left is immunohistochemistry of paraffin-embedded Human thyroid cancer tissue using CSB-PA199353(TNFRSF12A Antibody) at dilution 1/10, on the right is treated with fusion protein. (Original magnification: x200)

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.