





DNMT3L Antibody

Product Code	CSB-PA189054
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9UJW3
Immunogen	Fusion protein of Human DNMT3L
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA,IHC;ELISA:1:1000-1:2000,IHC:1:10-1:50
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	-20°C, pH7.4 PBS, 0.05% NaN ₃ , 40% Glycerol
Purification Method	Antigen affinity purification
Isotype	IgG
Alias	DNA (cytosine-5-)-methyltransferase 3-like
Immunogen Species	Homo sapiens (Human)
Target Names	DNMT3L

Image



The image on the left is immunohistochemistry of paraffin-embedded Human colon cancer tissue using CSB-PA189054(DNMT3L Antibody) at dilution 1/10, on the right is treated with fusion protein. (Original magnification: x200)



The image on the left is immunohistochemistry of paraffin-embedded Human thyroid cancer tissue using CSB-PA189054(DNMT3L Antibody) at dilution 1/10, on the right is treated with fusion protein. (Original magnification: x200)