




DOK1 Antibody

Product Code	CSB-PA040823
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q99704
Immunogen	Fusion protein of Human DOK1
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	ELISA,WB;ELISA:1:2000-1:5000,WB:1:500-1:2000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	-20°C, pH7.4 PBS, 0.05% NaN ₃ , 40% Glycerol
Purification Method	Antigen affinity purification
Isotype	IgG
Alias	docking protein 1, 62kDa (downstream of tyrosine kinase 1)
Immunogen Species	Homo sapiens (Human)
Target Names	DOK1

Image



Gel: 8%SDS-PAGE,Lysate: 40 µg,Lane 1-2:
Human hepatocellular carcinoma tissue, Jurkat
cells,Primary antibody: CSB-PA040823(DOK1
Antibody) at dilution 1/450 dilution,Secondary
antibody: Goat anti rabbit IgG at 1/8000
dilution,Exposure time: 20 seconds

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.