





TNFRSF25 Antibody

Product Code	CSB-PA213899
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q93038
Immunogen	Fusion protein of Human TNFRSF25
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA,IHC;ELISA:1:1000-1:2000,IHC:1:10-1:50
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	-20°C, pH7.4 PBS, 0.05% NaN ₃ , 40% Glycerol
Purification Method	Antigen affinity purification
Isotype	IgG
Alias	tumor necrosis factor receptor superfamily, member 25
Immunogen Species	Homo sapiens (Human)
Target Names	TNFRSF25

Image



The image on the left is immunohistochemistry of paraffin-embedded Human liver cancer tissue using CSB-PA213899(TNFRSF25 Antibody) at dilution 1/20, on the right is treated with fusion protein. (Original magnification: x200)



The image on the left is immunohistochemistry of paraffin-embedded Human cervical cancer tissue using CSB-PA213899(TNFRSF25 Antibody) at dilution 1/20, on the right is treated with fusion protein. (Original magnification: x200)

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.