




ELN Antibody

Product Code	CSB-PA591479
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P15502
Immunogen	Fusion protein of Human ELN
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, WB; ELISA:1:2000-1:5000, WB:1:400-1:2000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	-20°C, pH7.4 PBS, 0.05% NaN ₃ , 40% Glycerol
Purification Method	Antigen affinity purification
Isotype	IgG
Alias	elastin
Immunogen Species	Homo sapiens (Human)
Target Names	ELN

Image



Gel: 10%SDS-PAGE, Lysate: 40 µg, 60 µg,
Lane: Human hepatocellular carcinoma tissue,
Primary antibody: CSB-PA591479(ELN
Antibody) at dilution 1/400, Secondary antibody:
Goat anti rabbit IgG at 1/8000 dilution, Exposure
time: 2 minutes

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.