





ENPP5 Antibody

Product Code	CSB-PA970254
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9UJA9
Immunogen	Fusion protein of Human ENPP5
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	ELISA,IHC;ELISA:1:2000-1:10000,IHC:1:100-1:300
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	-20°C, pH7.4 PBS, 0.05% NaN ₃ , 40% Glycerol
Purification Method	Antigen affinity purification
Isotype	IgG
Alias	ectonucleotide pyrophosphatase/phosphodiesterase 5 (putative)
Immunogen Species	Homo sapiens (Human)
Target Names	ENPP5

Image



The image on the left is immunohistochemistry of paraffin-embedded Human breast cancer tissue using CSB-PA970254(ENPP5 Antibody) at dilution 1/60, on the right is treated with fusion protein. (Original magnification: x200)



The image on the left is immunohistochemistry of paraffin-embedded Human ovarian cancer tissue using CSB-PA970254(ENPP5 Antibody) at dilution 1/60, on the right is treated with fusion protein. (Original magnification: x200)

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.