





# DUSP14 Antibody

|                            |                                                               |
|----------------------------|---------------------------------------------------------------|
| <b>Product Code</b>        | CSB-PA055199                                                  |
| <b>Storage</b>             | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze. |
| <b>Uniprot No.</b>         | O95147                                                        |
| <b>Immunogen</b>           | Fusion protein of Human DUSP14                                |
| <b>Raised In</b>           | Rabbit                                                        |
| <b>Species Reactivity</b>  | Human,Mouse                                                   |
| <b>Tested Applications</b> | ELISA,IHC;ELISA:1:1000-1:2000,IHC:1:15-1:50                   |
| <b>Form</b>                | Liquid                                                        |
| <b>Conjugate</b>           | Non-conjugated                                                |
| <b>Storage Buffer</b>      | -20°C, pH7.4 PBS, 0.05% NaN <sub>3</sub> , 40% Glycerol       |
| <b>Purification Method</b> | Antigen affinity purification                                 |
| <b>Isotype</b>             | IgG                                                           |
| <b>Alias</b>               | dual specificity phosphatase 14                               |
| <b>Immunogen Species</b>   | Homo sapiens (Human)                                          |
| <b>Target Names</b>        | DUSP14                                                        |

## Image



The image on the left is immunohistochemistry of paraffin-embedded Human breast cancer tissue using CSB-PA055199(DUSP14 Antibody) at dilution 1/20, on the right is treated with fusion protein. (Original magnification: ×200)



The image on the left is immunohistochemistry of paraffin-embedded Human esophagus cancer tissue using CSB-PA055199(DUSP14 Antibody) at dilution 1/20, on the right is treated with fusion protein. (Original magnification: ×200)

## Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.